Gebruiksaanweisung

Schwenkriegel-System src

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Instructions for use

Swivel-type lock system src

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Thanks to prefabricated ceramic parts, the swivel-type lock system is especially suitable for the fabrication of dental restorations using CrCo alloys. The accurate impression of the ceramic parts facilitates finishing of the latch box after casting. The spring in the latch box ensures safe locking of the swivel-type lock when it is opened or closed.

2. Indication range

The swivel-type lock system src is used as a non-friction fixation element for removable dentures. It can also be used for monoreductors since the lock ensures reliable fixation. It is particularly suitable for the use in the one-piece casting technique.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Processing

The primary construction is waxed up in the usual way. The lock anchor is waxed to the primary construction at a previously determined angle using the paralleling mandrel and adapted towards the basal direction in accordance with the respective situation. The use of an oral transverse load distributor at the abutment crown is not required since the lock anchor features two grooves facing each other, which replace the distributor. Accordingly, enhanced aesthetics can be achieved when designing the crowns. When preparing bar constructions, the grooves can be filled with wax. After casting, the primary construction is fitted on the model. Lock anchor and the swivel-type lock are polished to high luster without any considerable dimensional changes. The ceramic lock blade must fit into the lock opening. The model is prepared for duplicating according to standard criteria. The primary construction must be blocked out to the basal direction. The opening of the lock is blocked out in a way to ensure that a shoulder of approx. 1 mm can be seen to ensure accurate positioning of the ceramic lock blade. Then the model is duplicated, e.g. using Exaktosil N 15 (REF 540 0103 4). The duplicate mould is filled with impression material. The ceramic lock blade is processed using a diamond grinding tool to ensure that it fits perfectly into the opening of the lock. Fit the ceramic axle in the ceramic lock blade and fix it; then both elements are integrated into the model (wax-up). The wax coating around the ceramic lock blade ensures that any reinforcement of this area is unnecessary. After casting, the ceramic elements are sandblasted using abrasive beads (50 µ) and the casting is sandblasted with aluminum oxide (50 µ) and fitted on. The opening for the steel pin is re-drilled using the Diatit Multidriill 1.0 (REF 330 0062 0) and milling and drilling oil (REF 550 0000 8) and the lock blade is inserted. Insert the spring with the open end first into the opening and place to the distal direction. Insert the lock blade and fix with the lock axle (see illustration). The lock blade is either fixed with rivets or resin to allow exchange of the spring at a later date.

4.1 Lock anchor

The upper jaw is the basis for proper placement of the lock anchors. Accordingly, the right lock anchor is placed in the right side of the upper jaw (1. quadrant) and the left lock anchor in the left side of the upper jaw. In the lower jaw, the right lock anchor is placed in the third quadrant and the left lock anchor in the fourth quadrant (mirror-inverted).

5. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

ad 4. Processing

Sectional view of the swivel-type lock system src

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